

Tork Blue Couch Roll Advanced



Article	152250
System	C1 - Couch roll system
Roll length	54.5 m
Roll width	50 cm
Roll diameter	12 cm
Number of sheets	165
Number of sheets	165
Sheet length	33 cm
Core inner diameter	5.1 cm
Ply	2
Embossing	Yes
Color	Blue

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in the healthcare environment. Rolls can be used with the Tork Couch Roll Dispenser.

www.tork.co.uk

Perforation for easy tear off and consumption control

Strong and easy to roll, for optimal protection

Better resistance and higher comfort due to 2-ply

Shipping data

	Consumer unit	Transport unit	Pallet
EAN	7310791230375	7310791230320	7322540716146
Pieces	1	9	162
Consumer units	-	9	162
Width	500 mm	500 mm	1161 mm
Height	120 mm	360 mm	1000 mm
Length	120 mm	360 mm	1200 mm
Volume	7.2 dm ³	64.8 dm ³	1.17 m ³
Net weight	980 g	8.8 kg	158.8 kg
Gross weight	1025 g	9.3 kg	167.5 kg
Packing material	Plastic	Plastic	-

Choose a dispenser**129184****Product certificates****Contact**

Essity UK Ltd
 Southfields Road
 Dunstable
 Bedfordshire
 LU6 3EJ

TorkCS.UK@essity.com

Tel 01582 677570

Fax 01582 677580